Virusanalytik

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ARNT EBINGER, OLAF DEGEN
BIOTECON DIAGNOSTICS GMBH, POTSDAM

Noroviruses (NoV) are the leading cause of non-bacterial gastroenteritis worldwide. Beside person-to-person transmission, food is considered as an important source of this viral infection. A novel Norovirus (GI, GII) plus Hepatitis A Virus Detection Kit (rRT-PCR) allows the simultaneous amplification and differentiation of NoV GI and II, hepatitis A virus and a process control virus in one reaction.

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
Literatur

[1] Kroneman A, Vega E, Vennema H et al. (2013) Proposal for a unified norovirus nomenclature and genotyping. *Arch Virol* 158:2059–2068

- [2] de Graaf M, van Beek J, Vennema H et al. (2015) Emergence of a novel GII.17 norovirus – End of the GII.4 area? *Euro Surveill* 20 (<http://dx.doi.org/10.2807/1560-7917.ES2015.20.26.21178>)
- [3] Centers for Disease Control and Prevention, Norovirus Worldwide – Global Trends, <https://www.cdc.gov/norovirus/worldwide.html>
- [4] Mäde D, Trübner K, Neubert E et al. (2013) Detection and typing of norovirus from frozen strawberries involved in a large-scale gastroenteritis outbreak in Germany. *Food Environ Virol*, doi: 10.1007/s12560-013-9118-0
- [5] Verhoef L, Hewitt J, Barclay L et al. (2015) Norovirus genotype profiles associated with foodborne transmission, 1999–2012. *Emerg Infect Dis* 21:592–599
- [6] International Organization for Standardization (2017) ISO/TS 15216-1:2017. Microbiology of the food chain – Horizontal method for determination of hepatitis A virus and norovirus using real-time RT-PCR, <https://www.iso.org/standard/65681.html>
- [7] Ebinger A (2017) foodproof Norovirus (GI, GII) plus Hepatitis A Virus Detection Kit, R 302 50. Validation Data Report, <http://www.bc-diagnostics.com>
- [8] Teunis PF, Moe CL, Liu P et al. (2008) Norwalk virus: how infectious is it? *J Med Virol* 80:1468–1476
- [9] Venter JM, van Heerden J, Vivier JC et al. (2007) Hepatitis A virus in surface water in South Africa: what are the risks? *J Water Health* 5:229–240
- [10] Bartsch C, Szabo K, Dinh-Thanh M et al. (2016) Comparison and optimization of detection methods for noroviruses in frozen strawberries containing different amounts of RT-PCR inhibitors. *Food Microbiol* 60:124–130
- [11] Scherer K, Mäde D, Ellerbroek L et al. (2009) Application of a Swab Sampling Method for the Detection of Norovirus and Rotavirus on Artificially Contaminated Food and Environmental Surfaces. *Food Environ Virol* 1:42

Korrespondenzadresse:

Arnt Ebinger
BIOTECON Diagnostics GmbH
Hermannswerder 17
D-14473 Potsdam
Tel.: 0331-2300-235
Fax: 0331-2300-299
aebinger@bc-diagnostics.com
www.bc-diagnostics.com

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© Springer-Verlag 2017